

# Bisulfite pyrosequencing protocol for Human sperm DNA

## Version 2

Zhaoxu Lu, Mei Qiang

### Abstract

**Citation:** Zhaoxu Lu, Mei Qiang Bisulfite pyrosequencing protocol for Human sperm DNA. **protocols.io**

dx.doi.org/10.17504/protocols.io.n52dg8e

**Published:** 30 Mar 2018

## Protocol

### Human Sperm DNA Isolation procedure

#### Step 1.

Extraction buffer of sperm DNA

AMOUNT

21.2 ml Additional info: 6mol/L guanidine thiocyanate

AMOUNT

600 µl Additional info: 5mol/L NaCl

AMOUNT

1 ml Additional info: 30% N-lauroylsarcosine sodium salt

AMOUNT

3 ml Additional info: 1mol/l dithiothreitol (DTT)

AMOUNT

600 µl Additional info: 10mg/ml proteinase K

AMOUNT

3.6 ml Additional info: Doubly deionized water

### Human Sperm DNA Isolation procedure

#### Step 2.

Add 150 µl of semen to 1.5 ml micro-centrifuge tube.

### Human Sperm DNA Isolation procedure

#### Step 3.

Wash with 1 ml of PBS (0.1 mol/L).

### Human Sperm DNA Isolation procedure

#### Step 4.

Centrifuge at 1500 ×g for 10 min at 4°C.

### Human Sperm DNA Isolation procedure

#### Step 5.

Repeat 2 times washing as step 3-4.

### Human Sperm DNA Isolation procedure

#### Step 6.

Add 0.5 ml extraction buffer into sperm pelle

### Human Sperm DNA Isolation procedure

**Step 7.**

Place in a 65°C water bath for 12 h.

Human Sperm DNA Isolation procedure

**Step 8.**

Cool at room temperature.

Human Sperm DNA Isolation procedure

**Step 9.**

Add 10 µl of RNase A (10 mg/ml), mix by pulse-vortexing for 15s, and incubate for 10 min at room temperature.

Human Sperm DNA Isolation procedure

**Step 10.**

Briefly centrifuge the tube.

Human Sperm DNA Isolation procedure

**Step 11.**

Add 510 µl of isopropanol and centrifuge at 10000 ×g for 10 min at 4°C.

Human Sperm DNA Isolation procedure

**Step 12.**

Add 800 µl of ethanol (75%), and reverse mixing for dozens of times.

Human Sperm DNA Isolation procedure

**Step 13.**

Incubate for 12 h at - 20°C.

Human Sperm DNA Isolation procedure

**Step 14.**

Centrifuge at 10000 ×g for 10 min at 4°C. Then dry sample at room temperature.

Human Sperm DNA Isolation procedure

**Step 15.**

Sperm DNA is dissolved in 50 µl of Elution Buffer.

Human Sperm DNA Isolation procedure

**Step 16.**

Incubate in a 65°C water bath for 2 h.

Procedure for bisulfite treatment

**Step 17.**

Add 130 µl of the CT Conversion Reagent solution to 1000 ng of your DNA sample in a PCR tube.

Procedure for bisulfite treatment

**Step 18.**

Place the sample tube in a thermal cycler and perform the following steps\*: ① 98°C for 10 minutes; ② 64°C for 2.5 hours; ③ 4°C.

Procedure for bisulfite treatment

**Step 19.**

Add 600 µl of M-Binding Buffer into a Zymo-Spin IC™ Column and place the column into a provided Collection Tube.

Procedure for bisulfite treatment

**Step 20.**

Load sample (from Step 2) into the Zymo-Spin IC™ Column containing the M-Binding Buffer. Close the cap and mix by inverting the column several times.

## Procedure for bisulfite treatment

### Step 21.

Centrifuge at full speed (>10,000 xg) for 30 seconds. Discard the flow-through.

## Procedure for bisulfite treatment

### Step 22.

Add 100 µl of M-Wash Buffer to the column. Spin at full speed for 30 seconds.

## Procedure for bisulfite treatment

### Step 23.

Add 200 µl of M-Desulphonation Buffer to the column and let stand at room temperature (20 °C - 30 °C) for 15-20 minutes. After the incubation, spin at full speed for 30 seconds.

## Procedure for bisulfite treatment

### Step 24.

Add 200 µl of M-Wash Buffer to the column. Spin at full speed for 30 seconds. Add another 200 µl of M-Wash Buffer and spin at top speed for 30 seconds.

## Procedure for bisulfite treatment

### Step 25.

Add 8 µl of M-Elution Buffer directly to the column matrix. Place the column into a 1.5 ml tube. Spin briefly at full speed to elute the DNA. Add 7 µl of M-Elution Buffer and additional repeated 1 time eluting was subsequently performed.

## Procedure for bisulfite treatment

### Step 26.

The DNA is ready for immediate analysis or can be stored at or below -20°C for later use.

## PCR amplification of bisulfite-treated sperm DNAs

### Step 27.

All reactions are performed with provided PCR mixtures (total volume at 25µl) provided in Table 1. Each reaction also contain 2.5 ul of CoralLoad Concentrate (10x) for checking amplicons on an agarose gel.

cmd [COMMAND \(Table 1 Components of PCR mixtures\)](#)

Components Volume (µl)	Final concentration	PyroMark PCR Master Mix, 2×	12.5	1×	CoraLoad Concentrate, 10×	2.5	1×	Q-solution, 5×	5	1×	primer forward (10 uM)	0.5	0.2uM	primer reverse (10 uM)	0.5	0.2uM	Template DNA	50ng	Final volume	25
------------------------	---------------------	-----------------------------	------	----	---------------------------	-----	----	----------------	---	----	------------------------	-----	-------	------------------------	-----	-------	--------------	------	--------------	----

## PCR amplification of bisulfite-treated sperm DNAs

### Step 28.

PCR and pyrosequencing primers are designed and listed in Table 2. Reverse primer is conjugated to biotin.

cmd [COMMAND \(Table 2 Primer sequences for PCR amplification\)](#)

DMR Forward primer	*Reverse primer	Sequencing primer	H19	GTATATGGGTATTTTTTGGAGGT
ATATCCTATTCCCAAATAA	TGGTTGTAGTTGTGGAAT	MEG3	GGGATTTTTGTTTTTTTTGTAGTAGG	
CCAACAAAACCCACCTATAAC	TTTGGGGTTGGGGTT	PEG3	TAATGAAAGTTTTGAGATTTGTTG	
CCTATAACAACCCACACCTATAC	GGGGGTAGTTGAGGTT			

## PCR amplification of bisulfite-treated sperm DNAs

### Step 29.

The PCR conditions are used as following: 94 °C for 15 min, followed by 45 cycles of 94 °C, 30 s, 56 °C, 30 s, 72°C, 30 s, and by a 72 °C final extension step for 10min.

## Pyrosequencing

### Step 30.

Add 40 µl of Binding Buffer, 3 µl of streptavidin-sepharose beads and 17 µl DDW into 20 µl of PCR

products.

#### Pyrosequencing

##### **Step 31.**

Seal film and shake at 1400 rpm for 10 min at room temperature.

#### Pyrosequencing

##### **Step 32.**

PCR products on streptavidin-sepharose beads are washed with ethanol (10%) for 5s.

#### Pyrosequencing

##### **Step 33.**

Place sample (from step 3) into denaturation solution for 5s.

#### Pyrosequencing

##### **Step 34.**

Place sample (from step 4) into Wash Buffer for 10s for getting purified biotinylated single stranded PCR products. These single stranded PCR products are isolated using the Pyrosequencing Work Station.

#### Pyrosequencing

##### **Step 35.**

Transfer purified biotinylated single stranded PCR products into PSQ 96 Plate Low with 40  $\mu$ l annealing buffer and 1.6  $\mu$ l sequencing primer (10 $\mu$ mol/L).

#### Pyrosequencing

##### **Step 36.**

Heat PSQ 96 Plate Low at 80 °C for 2 minutes.

#### Pyrosequencing

##### **Step 37.**

Undergo pyrosequencing on a Pyromark Q96 MD pyrosequencing instrument and sequence using PyroMark Gold Q96 kit.

#### Pyrosequencing

##### **Step 38.**

The degree of methylation at each CpG site is determined using PyroMark CpG Software (Biotage AB, Uppsala, Sweden). Pyrosequencing assays are performed in duplicate in sequential runs (technical replicates), and the values show represent the mean methylation for each individual CpG site.